



**FRESENIUS  
KABI**



# Aurora Plasmapheresis System

Experience the value of donor-focused design.

From the moment donors arrive to the moment they leave, every detail makes a difference in the donor experience. How donors feel when they're connected to a device. How smoothly and safely the process goes. How confidently your phlebotomists and technicians respond. It's what keeps donors coming back. And it's all been considered in the Aurora system.

---

Optimized Nomogram collects more plasma per donation

---

More productive procedures mean that your teams and donors can help more patients. So on Aurora we've added an Optimized Nomogram feature to increase collection yield. With automated calculation of maximum plasma volume based on weight and hematocrit, you'll never miss out on results that would be more rewarding, for everyone.

# Aurora Plasmapheresis System

The foundation of the Aurora system is our proprietary Spinning Membrane technology – a rapidly rotating membrane with tangential filtration that enables the collection of low cell count plasma.

## Let them donate in comfort

Whether first-time or long-time donors, they want every procedure to go just right.

**Convenient side displays** let donors monitor their own progress.

**Adaptive Vein Control** monitors venous pressure and automatically adjusts draw and return rates.

**Less than 200 mL ECV** – the lowest in the industry – helps reduce the potential of future deferrals.



## Free to focus on donors


When your phlebotomists are confident in the tools they use, they're able to give donors an extra dose of their attention.

**Intuitive touch screen interfaces and visual prompts** support smooth operation.

**SpikeSmart® system** promotes safety by reducing the risk of interchanging saline and AC connections.

**Instrument Signal Light** enables operators to see changes in the donor procedure from a distance.

**Wifi-enabled devices** allow for bi-directional communication with the DXT® Data Management System.



## Every device at the ready

You never want to make a donor wait, or turn them away, because you don't have a device ready to go.

**Preventative maintenance routines** are quick and easy, taking just minutes to keep devices up and running.

**Easy-access back panel** facilitates quick part repairs and replacements.

**Customer satisfaction** comes first. Our plasma consultants and service and support teams are always ready to help your center.



Contact your local Fresenius Kabi representative to learn more about the Aurora Plasmapheresis system.

This marketing reflects compliance with the applicable CE Marketing requirements for medical devices.



Refer to the Operator's Manuals for a complete list of warnings and cautions associated with the use of this product. All trademarks shown are the property of their respective owners.



**FRESENIUS  
KABI**

Fresenius Kabi AG  
Else-Kröner-Straße 1  
61352 Bad Homburg  
Germany

Phone: +49 61 72 608-0

[www.fresenius-kabi.com](http://www.fresenius-kabi.com)